


ACCENT: Transferring Vocabulary Files to & From the Device & NuVoice PASS Software



 Click [HERE](#) to download NuVoice PASS

Technical Support:

Tel: +44 (0)1733 370 470 (Option 2)

Email: info@liberator.co.uk

Regional AAC Consultants:

If you wish to speak to your Regional AAC Consultant, please call +44 (0)1733 370 470 (Option 3)

This quick reference guide will show you how to save and import an individual vocabulary or ALL vocabularies to and from a device, and the NuVoice PASS software.

Saving an Individual Vocabulary

1. Insert the memory stick into the device/PC
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
3. Select **IMPORT/SAVE... > SAVE A VOCABULARY**.

4. Select the user area to save (this will be one of the 6 listed, usually marked **CURRENT AREA**).
5. Choose the location of your save by selecting **CHOOSE ANOTHER FOLDER** and selecting your memory stick as your location, then select **OK**.
6. Type the name for your save (e.g. Client NameDDMMYY).
7. Select **OK** and the **WAIT (Transfer in Progress)** dialogue box appears.
8. Once complete, select **OK** several times to return to the home screen.

Importing an Individual Vocabulary

1. Insert the memory stick into the device/PC.
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
3. Select **IMPORT/SAVE... > IMPORT A VOCABULARY**.
4. Select the area in which you want to load the new vocabulary into.
5. Choose the location of your file by using **CHOOSE ANOTHER FOLDER** and selecting your memory stick as the source.
6. Select the file to load and the **WAIT (Transfer in Progress)** dialogue box appears.
7. Once complete, select **OK** several times to return to the home screen.

Saving ALL Vocabularies

1. Insert the memory stick into the device/PC.
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
3. Select **IMPORT/SAVE... > SAVE ALL VOCABULARIES**.
4. Choose the location of your save by selecting **CHOOSE ANOTHER FOLDER** and selecting your memory stick as your location, then select **OK**.
5. Type the name for your save (e.g. Client NameDDMMYY).
6. Select **OK** and the **WAIT (Transfer in Progress)** dialogue box appears.
7. Once complete, select **OK** several times to return to the home screen.

Importing ALL Vocabularies

1. Insert the memory stick into the device/PC.
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).

3. Select **IMPORT/SAVE... > IMPORT ALL VOCABULARIES > YES**
4. Choose the location of your file by using **CHOOSE ANOTHER FOLDER** and selecting your memory stick as the source.
5. Select the file to load and the **WAIT (Transfer in Progress)** dialogue box appears.
6. Once complete, select **OK** several times to return to the home screen.

Switching to the Correct Vocabulary

1. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
2. Select **VOCABULARY OPTIONS > CHANGE VOCABULARY**.
3. Choose the vocabulary you want to switch to
4. Press **OK > GOTO HOME**

Saving a Page

1. Insert the memory stick into the device/PC.
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
3. Select **IMPORT/SAVE... > SAVE PAGES & ACTIVITIES**.
4. Select **SPELL PAGE TO FIND**.
5. Type the name of your page and select **OK**.
6. Choose the page to save.
7. Choose the location of your save by selecting **CHOOSE ANOTHER FOLDER** and selecting your memory stick as your location, then select **OK**.
8. Select **OK** to save.

Importing a Page

1. Insert the memory stick into the device/PC.
2. Select **TOOLS > TOOLBOX** (or press F9 on your keyboard if using the Emulation Software).
3. Select **IMPORT/SAVE... > IMPORT PAGES & ACTIVITIES**.
4. Choose the location of your save by selecting **CHOOSE ANOTHER FOLDER** and selecting your memory stick as your location, then select **OK**.
5. Wait, and once complete select **OK** several times to return to the home screen.