


ACCENT: Using LAMP Words For Life



 Click [HERE](#) to download NuVoice PASS

Technical Support:

Tel: +44 (0)1733 370 470 (Option 2)

Email: info@liberator.co.uk

Regional AAC Consultants:

If you wish to speak to your Regional AAC Consultant, please call +44 (0)1733 370 470 (Option 3)

Loading LAMP Words For Life Vocabularies

1. Select **Tools > Toolbox > Vocabulary Options**.
2. Select **Replace Vocabulary** and select the area to replace.
NOTE: replacing a user area will **OVERWRITE** the existing area. Be sure the area is saved if there is content you don't want to lose!
3. Select **LAMP Words for Life**.
4. Select **1-Hit, Transition, or Full** and select **YES**.

5. Type **replace** and select **OK**. Select **YES** when prompted to switch to this user area.
6. Select **OK**, then **GOTO HOME** to exit.

Changing the Voice

1. Select **Tools > Toolbox > Speech**.
2. Select **Change Voice**.
3. Select the desired voice and wait for it to load. To listen to the voice, select **Test Speech**.
4. Adjust **rate, volume, bass, and treble** as needed.

Pronunciation Dictionary

"Teach" the device to pronounce a word correctly.

1. Select **Tools > Toolbox > Pronunciations**.
2. Choose **Main Voice**.
3. Select **Add Pronunciation**.
4. Enter the correct spelling of the word. Select **OK**.
5. Enter the "**phonetic spelling**" of the word. To test the pronunciation, select the word added and press **Test Pronunciation**. Select **OK**, then **GOTO HOME** to exit.

Vocabulary Backup

1. Insert a USB flash drive into the USB port.
2. Select **Tools > Toolbox > Import/Save**.
3. Select **Save a Vocabulary** followed by the area you want to save, or **Save ALL Vocabularies** to save the global NuVoice settings, including IR.
4. Make sure the pre-set folder is the USB drive and select **OK**. If not, select **Choose Another Folder** and navigate to the USB drive. Select **OK** and type the name for the file.
5. Select **OK** and wait for the transfer. Select **OK** when finished.

Disabling the Toolbox

To disable the toolbox:

1. Select **Tools > Toolbox > Advanced Settings**.
2. Select **System Lock Settings**.
3. Turn **Toolbox Restrictions** ON and DISABLE the **On-Screen Tools Key**.

4. Select whether restriction mode is **USB LOCK** or **PASSWORD LOCK**.
 - a. If **USB LOCK**, insert a USB flash drive into the USB port and select **Create USB Lock Override**.
 - b. If **PASSWORD LOCK**, select **Create Toolbox Password**, enter a password and select **OK** to disable the toolbox. *Remember your password!*

To enable the toolbox:

If locked with a password:

1. Press the **Tools** button at the bottom of the device.
2. Enter your **password** and select **OK**.
3. Select **Toolbox > Advanced Settings**.
4. Select **System Lock Settings** and DISABLE the **Toolbox Restrictions**.

TIP: If you forget your password the toolbox can be unlocked with a USB flash driving containing the USB Lock Override file.

If locked with USB lock:

1. Insert USB flash drive containing the USB Lock Override file.
2. Press the **Tools** button at the bottom of the device.
3. Follow above instructions from step 3.

Additional Options:

Within the System Lock Settings menu, you can Enable or Disable the following buttons:

- Hardware Volume Keys, Tools Key (not recommended), On-Screen Tools Key, Help Button, Pulldown Menu Button, and the Left, Right & Bottom On-Screen Tools Buttons.

Taking a Photo

1. Select **Tools > Toolbox > Camera**.
2. Frame the picture and select **Take Picture**.
3. Select **Save as Icon**. Rename the picture if you want.
4. Select **OK** followed by **Cancel**.
Note: The above steps will place the picture in the Imports category when you select **Change Picture** (See step 2b under **Setup Key**).

Setup Key

1. Be sure you're looking at the button you want to edit. Select **Tools**, then **Setup Key**. Touch the button you wish to edit.
2. Select **Change Picture** and either:
 - a. **Spell Icon to Find**: Spell the beginning letters for the icon and select OK. Touch the icon you want.
 - b. Select by **Category**: Look through the categories for an icon. **NOTE: Camera pictures are in Imports.**
3. Look to see if "**Message to Speak...**" is correct. If not, select **Add Message or Insert Tool** to type text that will appear in the text area and be spoken. Be sure to put a space after the last character. Select **OK**.
4. Check to see if the label is correct. If not, select **Change Label**, type the new label, and select **OK**.
5. Select **OK** to exit, or choose **Edit Another Button**.

Creating and Linking Pages

1. Be sure you're looking at the button you want to edit. Select **Tools** > **Setup Key**. Touch the button you wish to link to a new page.
2. Select **Change Picture** to search for a picture for the button.
3. Select **Add Message or Insert Tool**.
4. Select **Clear** followed by **Page Link** and **CLOSE**.
5. If the page is already created, find and select it from the menu. If it is not already created, select **Create New Page**, type the name for the page, select **OK**, and choose keyboard size 84. Select **OK** two times.
6. Select **Change Label** to change the label on the key. Labels for page links are typically in all CAPS. Select **OK**.
7. Select **OK** to exit when finished.

Vocabulary Builder

Vocabulary Builder is a therapeutic tool designed to teach new words within the LAMP Words for Life® vocabulary. The strategy of "masking" vocabulary allows a small amount of vocabulary to be taught on the device while keeping learned motor plans constant.


Turn Vocabulary Builder On/Off

1. Select the **Tools** button.
2. Select **Vocabulary Builder On/Off**. When Vocabulary Builder is ON, a vocabulary overlay will appear in the status display.



How to create a set of words:

1. Select the **Tools** button.
 2. Select **Vocabulary Builder Menu**.
 3. Select **Mask All***.
 4. Select **Spell to Unmask**.
 5. Type the word(s) to unmask using commas to separate words.
 6. Select **OK**, then,
 7. Select **OK** to exit the Vocabulary Builder Menu.
- * Skip this step if you are adding words to an existing set.



Saving Sets:

If you intend to keep the sequence sets you have built, you need to SAVE them!!

1. Select the **Tools** button.
2. Select **Vocabulary Builder Menu**.
3. Select **Save as Sequence Set**.
4. Type the name for the sequence set.
5. Select **OK** 2 times.

Loading Sets:

After creating a sequence set, it can be retrieved from the VOCABULARY BUILDER MENU.

1. Select the **Tools** button.
2. Select **Vocabulary Builder Menu**.
3. Select **Advanced Options**.
4. Select **Use Sequence Set**.
5. Select **Load Sequence Set**.
6. Choose the desired set.
7. Select **OK** 2 times.

Active Side When Spelling

When adding additional words to an existing set, make sure the right side has a tan background. Tap the right side to turn it from gray to tan.

Remember to follow the child's lead!

Use sequence sets to isolate words you intend to teach in structured or one-on-one activities. Add or change words frequently as the activity changes and turn Vocabulary Builder on/off as needed during unstructured time so the individual has an opportunity to independently explore. Remember to choose words that can easily generalise from one activity to another.